

Rapid Extraction Protocol

Contributed by Gina Brown-Guedira
USDA-ARS Genotyping Facility
Manhattan, KS

This protocol has been submitted for publication as a NOTE in Crop Science

1. Add coleoptiles or leaf tissue directly into a 96-well flat bottom tissue culture plate. Roll up tissue to fit into well, do not let tissue stick to the side of the plate. Freeze tissue at -20°C until ready to extract. **Fresh tissue works best.
2. Add 40ul 0.25 NaOH directly into wells containing leaf tissue.
3. Heat microtiter plate in 95°C water bath for 1 min.
4. Put stainless steel cylinders on top of tissue. Using the matrix mill, grind tissue for 7-10 min.
5. After grinding remove cylinders. Add 130ul of 0.1 Tris-HCl to each well.
6. Spin plate 3500rpm, 10min.
7. Remove 150ul supernatant and put into a 200uL 96-well conical bottom PCR plate
8. Add 15ul 3M NaOAc and 120ul 100% isopropanol
9. Freeze -80°C 1hr. or -20°C overnight
10. Centrifuge 3500rpm, 30min
11. Remove isopropanol and centrifuge upside down 600rpm, 1min
12. Add 200ul 70% EtOH
13. Centrifuge 3500rpm, 20min
14. Remove EtOH and centrifuge upside down 600rpm, 1min
15. Resuspend in 20ul TE buffer

Products we use:

Germination containers

8-well rectangular tissue culture dish (Nunc, Cat. No. 176600)

Extraction plastics

96-well flat bottom tissue culture dish (Nunc, Cat No. 167008)

Conical bottom 96-well PCR plate (Genemate®, Cat no. T-3031-21)

Mention of trademark or proprietary product does not constitute a guarantee or warranty by the USDA and does not imply its approval over other suitable products.

Questions? Contact Dr. Kristi Hill-Ambroz, e-mail: wheatgenes@yahoo.com

<http://maswheat.ucdavis.edu>